

# mattress

for behavioral health



# Mattress for Behavioral Health

A resilient mattress designed to enhance patient safety and comfort, featuring 5 distinct pressure reduction zones, cooling technology and inverted welded seams.

## Features

- Seams designed to resist bed bugs and contaminants
- 5 distinct pressure reduction zones to increase patient comfort
- Cooling technology helps dissipate heat
- Air flow through zoned channels improves air circulation

## Options

- Custom sizes available

## Construction

- TPU nylon fabric cover
- Inverted welded seams
- Base layer of foam is 3.5" thick and has 5 distinct pressure reduction zones
- Top layer is 1.5" high energy foam
- 3" wide extra firm side rails on each side to prevent edge collapse
- Required fire barrier is placed over the mattress core in such a way that it completely encapsulates the entire mattress forming a full fire barrier.
- Flame retardant (CFR 16 part1632 and 1633, BFD IX-11, TB 129, TB 603, TB121, NFPA 701, and TB117)
- 100% recyclable foam core
- Manufactured in the USA

## Cleaning Instructions

Use mild soap suds and lukewarm water, rinse well with water and allow to dry. **DO NOT USE HARSH CLEANSERS, SOLVENTS OR DETERGENTS.** For hard to clean spots use standard vinyl cleaners and soft bristle brush. To disinfect dilute disinfectants and/or germicides as specified on the product label.

Model number: SF36805

Length: 80"

Width: 36"

Height: 5"



Model number: SF36805-PW

Length: 80"

Width: 36"

Height: 5"



## Casegoods for behavioral health



Frontier



Fortress



Gibraltar



Liberty



STANCE | HEALTHCARE

45 Goodrich Drive Kitchener, ON N2C 0B8 tel 1.877.395.2623 fax 1.877.395.6270 stancehealthcare.com

S04-19 © Stance Healthcare Inc. all rights reserved.

